Basic techniques with Quick Band Professional

The training can also be carried without a quick band

<u>Contents</u>: 5 min jumping rope, 10 min gymnastics, 3 min specific warm-up, 40-45 min technique exercises including cool down. Total approx. 60 min.

Additional options: Individual stretching / warmup before the unit is recommended.

Required equipment: jump rope, Quick Band Professional (alternatively other bands), gym mat

Jump rope exercise, mat, Quick Band, which is attached to the upper body











Alternatives:

- Jump rope: stepping, covering u. To deceive
- Exercise mat: clean floor
- Quick-Band: Other bands (on the foot or body) or only techniques without an additional device

Note:

• If you train without an additional device, then train individually with more speed.

Rough process

5 min: jumping rope. Both legs, one-legged, cross turns, double jumps, etc ...

10 min: stretching

3-5 min: put on quickband and Present alternatives.

7 min: specific warm up 8 repetitions 20 sec break. 3 series (time 10/20 * 3)

Kneeling left / right 8 WH * 3 series (I 8 * 2)

Side step knee pull

Approach split right, set down, backstitch 8 WH * 2

Neck approach left, set down at front, 8 WH * 2 step

Cut-split approach, lower right, backstitch 8 WH * 2

Cut-split approach, set off left front, 8 WH * 2 step

2-3 min: break and Explanations
30 min: techniques with / without Quickband (10/20 * 2)
Push kick right / left (2 min)
Cut-Kick R / L (2 min)
Break (2 min)
Paltung-chagi R / L (2 min)
Cut-splitting R / L (2 min)
Break (3 min)
Faust, block in front R / L (2 min)
Faust, block back R / L (2 min)
Break (2 min)
Free kicking front leg R 8/20 * 2
Free kicking front leg L 8/20 * 2 (2 min)
Break (2 min)
Free kicking all techniques 8/20 * 4 (2 min)
Break (2 min)
Free kicking all techniques (2 min)
5 min: cool down with easy stitching and To deceive
Total approx. 55 - 65 min. With previous independent warm-up and If necessary, further free rounds can be carried out for about 1.5 hours.

Georg Streif

With sporty greetings